

Package ‘PKNCA’

June 13, 2018

Type Package

Title Perform Pharmacokinetic Non-Compartmental Analysis

Version 0.8.5

Imports dplyr (>= 0.5.0), digest, plyr, nlme, parallel, lattice,
rlang, stats, tidyr, utils

Suggests cowplot, ggplot2, knitr, rmarkdown, testthat

Description Compute standard Non-Compartmental Analysis (NCA) parameters and summarize them. In addition to this core work, it also provides standardized plotting routines, basic assessments for biocomparison or drug interaction, and model-based estimation routines for calculating doses to reach specific values of AUC or Cmax.

License AGPL-3

URL <https://github.com/billdenney/pknca>

BugReports <https://github.com/billdenney/pknca/issues>

NeedsCompilation no

VignetteBuilder knitr

RoxygenNote 6.0.1

Author Bill Denney [aut, cre],
Clare Buckeridge [aut],
Sridhar Duvvuri [ctb]

Maintainer Bill Denney <wdenney@humanpredictions.com>

Repository CRAN

Date/Publication 2018-06-13 07:55:49 UTC

R topics documented:

add.interval.col	4
addProvenance	6
adj.r.squared	6

AIC.list	7
as.data.frame.PKNCAresults	7
business.mean	8
check.conc.time	9
check.conversion	10
check.interval.deps	10
check.interval.specification	11
checkProvenance	11
choose.auc.intervals	12
clean.conc.blq	13
clean.conc.na	14
exclude	15
exclude_nca	16
find.tau	17
findOperator	18
fit_half_life	18
formula.parseFormula	19
formula.PKNCAconc	19
formulaOrNames	20
geomean	20
get.best.model	21
get.first.model	22
get.interval.cols	22
get.parameter.deps	23
getAttributeColumn	23
getColumnValueOrNot	24
getData.PKNCAconc	24
getData.PKNCAdata	25
getData.PKNCAresults	25
getDataName.PKNCAconc	26
getDepVar	26
getGroups.PKNCAconc	27
getIndepVar	28
interp.extrap.conc	28
merge.splitlist	30
model.frame.PKNCAconc	31
normalize_exclude	32
parseFormula	32
pk.business	33
pk.calc.ae	34
pk.calc.aucint	34
pk.calc.aucpext	36
pk.calc.auxc	36
pk.calc.c0	39
pk.calc.cav	40
pk.calc.ceoi	40
pk.calc.cl	41
pk.calc.clast.obs	42

pk.calc.clr	42
pk.calc.cmax	43
pk.calc.ctrough	43
pk.calc.deg.fluc	44
pk.calc.dn	44
pk.calc.f	45
pk.calc.fe	45
pk.calc.half.life	46
pk.calc.kel	47
pk.calc.mrt	48
pk.calc.mrt.md	49
pk.calc.ptr	49
pk.calc.swing	50
pk.calc.thalf.eff	50
pk.calc.tlag	51
pk.calc.tlast	51
pk.calc.tmax	52
pk.calc.vd	53
pk.calc.vss	53
pk.calc.vz	54
pk.nca	54
pk.nca.interval	55
pk.tss	56
pk.tss.data.prep	57
pk.tss.monoexponential	57
pk.tss.monoexponential.individual	58
pk.tss.monoexponential.population	59
pk.tss.stepwise.linear	60
PKNCA	61
PKNCA.choose.option	61
PKNCA.options	62
PKNCA.options.describe	63
PKNCA.set.summary	64
PKNCAconc	65
PKNCAdata	66
PKNCAdose	67
PKNCAresults	68
plot.PKNCAconc	69
print.PKNCAconc	70
print.PKNCAdata	70
print.provenance	71
roundingSummarize	71
roundString	72
setAttributeColumn	72
setDuration	73
setExcludeColumn	74
setRoute	74
signifString	75

sort.interval.cols	76
split.PKNCAconc	76
summary.PKNCAdata	77
summary.PKNCAresults	77
superposition	78
tss.monoexponential.generate.formula	80

Index	81
--------------	-----------

add.interval.col	<i>Add columns for calculations within PKNCA intervals</i>
------------------	------------------------------------------------------------

Description

Add columns for calculations within PKNCA intervals

Usage

```
add.interval.col(name, FUN, values = c(FALSE, TRUE), depends = c(),
  desc = "", formalsmap = list(), datatype = c("interval", "individual",
  "population"))
```

Arguments

name	The column name as a character string
FUN	The function to run (as a character string) or NA if the parameter is automatically calculated when calculating another parameter.
values	Valid values for the column
depends	Character vector of columns that must be run before this column.
desc	A human-readable description of the parameter (<=40 characters to comply with SDTM)
formalsmap	A named list mapping parameter names in the function call to NCA parameter names. See the details for information on use of formalsmap.
datatype	The type of data used for the calculation

Details

The formalsmap argument enables mapping some alternate formal argument names to parameters. It is used to generalize functions that may use multiple similar arguments (such as the variants of mean residence time). The names of the list should correspond to function formal parameter names and the values should be one of the following:

- For the current interval:
 - character strings of NCA parameter name** The value of the parameter calculated for the current interval.
 - "conc"** Concentration measurements for the current interval.

"time" Times associated with concentration measurements for the current interval (values start at 0 at the beginning of the current interval).

"volume" Volume associated with concentration measurements for the current interval (typically applies for excretion parameters like urine).

"duration.conc" Durations associated with concentration measurements for the current interval.

"dose" Dose amounts associated with the current interval.

"time.dose" Time of dose start associated with the current interval (values start at 0 at the beginning of the current interval).

"duration.dose" Duration of dose (typically infusion duration) for doses in the current interval.

"route" Route of dosing for the current interval.

"start" Time of interval start.

"end" Time of interval end.

"options" PKNCA.options governing calculations.

- For the current group:

"conc.group" Concentration measurements for the current group.

"time.group" Times associated with concentration measurements for the current group (values start at 0 at the beginning of the current interval).

"volume.group" Volume associated with concentration measurements for the current interval (typically applies for excretion parameters like urine).

"duration.conc.group" Durations associated with concentration measurements for the current group.

"dose.group" Dose amounts associated with the current group.

"time.dose.group" Time of dose start associated with the current group (values start at 0 at the beginning of the current interval).

"duration.dose.group" Duration of dose (typically infusion duration) for doses in the current group.

"route.group" Route of dosing for the current group.

Value

NULL (Calling this function has a side effect of changing the available intervals for calculations)

Examples

```
## Not run:
add.interval.col("cmax",
  FUN="pk.calc.cmax",
  values=c(FALSE, TRUE),
  desc="Maximum observed concentration",
  depends=c())
add.interval.col("cmax.dn",
  FUN="pk.calc.dn",
  values=c(FALSE, TRUE),
  desc="Maximum observed concentration, dose normalized",
  formalsmap=list(parameter="cmax"),
```

```

depends=c("cmax"))

## End(Not run)

```

addProvenance	<i>Add a hash and associated information to enable checking object provenance.</i>
---------------	------------------------------------------------------------------------------------

Description

Add a hash and associated information to enable checking object provenance.

Usage

```
addProvenance(object, replace = FALSE)
```

Arguments

object	The object to add provenance
replace	Replace provenance if the object already has a provenance attribute. (If the object already has provenance and replace is FALSE, then an error will be raised.)

Value

The object with provenance as an added item

See Also

[checkProvenance](#)

adj.r.squared	<i>Calculate the adjusted r-squared value</i>
---------------	-----------------------------------------------

Description

Calculate the adjusted r-squared value

Usage

```
adj.r.squared(r.sq, n)
```

Arguments

r.sq	The r-squared value
n	The number of points

Value

The numeric adjusted r-squared value

AIC.list	<i>Assess the AIC for all models in a list of models</i>
----------	----------------------------------------------------------

Description

Assess the AIC for all models in a list of models

Usage

```
## S3 method for class 'list'
AIC(object, ..., assess.best = TRUE)
```

Arguments

object	the list of models
...	parameters passed to the underlying AIC function (typically the parameter k)
assess.best	determine which model is the best (by lowest AIC)

Value

a data frame with row names matching the names of the list `x` and columns for degrees of freedom (df) and AIC. If `assess.best` is true, then there will be another column `isBest`.

See Also

[get.best.model](#)

as.data.frame.PKNCAresults	<i>Extract the parameter results from a PKNCAresults and return them as a data frame.</i>
----------------------------	-------------------------------------------------------------------------------------------

Description

Extract the parameter results from a PKNCAresults and return them as a data frame.

Usage

```
## S3 method for class 'PKNCAresults'
as.data.frame(x, ..., out.format = c("long", "wide"))
```

Arguments

x	The object to extract results from
...	Ignored (for compatibility with generic as.data.frame)
out.format	Should the output be 'long' (default) or 'wide'?

Value

A data frame of results

business.mean	<i>Generate functions to do the named function (e.g. mean) applying the business rules.</i>
---------------	---------------------------------------------------------------------------------------------

Description

Generate functions to do the named function (e.g. mean) applying the business rules.

Usage

```
business.mean(x, ...)
```

```
business.sd(x, ...)
```

```
business.cv(x, ...)
```

```
business.geomean(x, ...)
```

```
business.geocv(x, ...)
```

```
business.min(x, ...)
```

```
business.max(x, ...)
```

```
business.median(x, ...)
```

```
business.range(x, ...)
```

Arguments

x vector to be passed to the various functions

... Additional arguments to be passed to the underlying function.

Value

The value of the various functions or NA if too many values are missing

Functions

- `business.sd`: Compute the standard deviation with business rules.
- `business.cv`: Compute the coefficient of variation with business rules.
- `business.geomean`: Compute the geometric mean with business rules.
- `business.geocv`: Compute the geometric coefficient of variation with business rules.

- `business.min`: Compute the minimum with business rules.
- `business.max`: Compute the maximum with business rules.
- `business.median`: Compute the median with business rules.
- `business.range`: Compute the range with business rules.

See Also

`pk.business`

<code>check.conc.time</code>	<i>Verify that the concentration and time are valid</i>
------------------------------	---------------------------------------------------------

Description

If the concentrations or times are invalid, will provide an error. Reasons for being invalid are

- `time` is not a number
- `conc` is not a number
- Any `time` value is NA
- `time` is not monotonically increasing
- `conc` and `time` are not the same length

Usage

```
check.conc.time(conc, time, monotonic.time = TRUE)
```

Arguments

<code>conc</code>	Measured concentrations
<code>time</code>	Time of the measurement of the concentrations
<code>monotonic.time</code>	Must the time be unique and monotonically increasing?

Details

Some cases may generate warnings but allow the data to proceed.

- A negative concentration is often but not always an error; it will generate a warning.

Value

None

check.conversion	<i>Check that the conversion to a data type does not change the number of NA values</i>
------------------	-----------------------------------------------------------------------------------------

Description

Check that the conversion to a data type does not change the number of NA values

Usage

```
check.conversion(x, FUN, ...)
```

Arguments

x	the value to convert
FUN	the function to use for conversion
...	arguments passed to FUN

Value

FUN(x, ...) or an error if the set of NAs change.

check.interval.deps	<i>Take in a single row of an interval specification and return that row updated with any additional calculations that must be done to fulfill all dependencies.</i>
---------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------

Description

Take in a single row of an interval specification and return that row updated with any additional calculations that must be done to fulfill all dependencies.

Usage

```
check.interval.deps(x)
```

Arguments

x	A data frame with one or morw rows of the PKNCA interval
---	----------------------------------------------------------

Value

The interval specification with additional calculations added where requested outputs require them.

See Also

[check.interval.specification](#)

`check.interval.specification`*Check the formatting of a calculation interval specification data frame.*

Description

Calculation interval specifications are data frames defining what calculations will be required and summarized from all time intervals. Note: parameters which are not requested may be calculated if it is required for (or computed at the same time as) a requested parameter.

Usage

`check.interval.specification(x)`

Arguments

`x` The data frame specifying what to calculate during each time interval

Details

`start` and `end` time must always be given as columns, and the `start` must be before the `end`. Other columns define the parameters to be calculated and the groupings to apply the intervals to.

Value

`x` The potentially updated data frame with the interval calculation specification.

See Also

[check.interval.deps](#), [get.parameter.deps](#), [get.interval.cols](#)

`checkProvenance`*Check the hash of an object to confirm its provenance.*

Description

Check the hash of an object to confirm its provenance.

Usage

`checkProvenance(object)`

Arguments

`object` The object to check provenance for

Value

TRUE if the provenance is confirmed to be consistent, FALSE if the provenance is not consistent, or NA if provenance is not present.

See Also

[addProvenance](#)

choose.auc.intervals *Choose intervals to compute AUCs from time and dosing information*

Description

Intervals for AUC are selected by the following metrics:

1. If only one dose is administered, use the `PKNCA.options("single.dose.auc")`
2. If more than one dose is administered, estimate the AUC between any two doses that have PK taken at both of the dosing times and at least one time between the doses.
3. For the final dose of multiple doses, try to determine the dosing interval (τ) and estimate the AUC in that interval if multiple samples are taken in the interval.
4. If there are samples $> \tau$ after the last dose, calculate the half life after the last dose.

Usage

```
choose.auc.intervals(time.conc, time.dosing, options = list(),
  single.dose.aucs = NULL)
```

Arguments

<code>time.conc</code>	Time of concentration measurement
<code>time.dosing</code>	Time of dosing
<code>options</code>	List of changes to the default <code>PKNCA.options</code> for calculations.
<code>single.dose.aucs</code>	The AUC specification for single dosing.

Value

A data frame with columns for `start`, `end`, `auc.type`, and `half.life`. See [check.interval.specification](#) for column definitions. The data frame may have zero rows if no intervals could be found.

See Also

[pk.calc.auc](#), [pk.calc.aumc](#), [pk.calc.half.life](#), [find.tau](#), [check.interval.specification](#), [PKNCA.options](#)

clean.conc.blq	<i>Handle BLQ values in the concentration measurements as requested by the user.</i>
----------------	--------------------------------------------------------------------------------------

Description

Handle BLQ values in the concentration measurements as requested by the user.

Usage

```
clean.conc.blq(conc, time, ..., options = list(), conc.blq = NULL,
               conc.na = NULL, check = TRUE)
```

Arguments

conc	Measured concentrations
time	Time of the concentration measurement
...	Additional arguments passed to <code>clean.conc.na</code>
options	List of changes to the default PKNCA.options for calculations.
conc.blq	How to handle a BLQ value that is between above LOQ values? See details for description.
conc.na	How to handle NA concentrations. (See clean.conc.na)
check	Run check.conc.time ?

Details

NA concentrations (and their associated times) will be handled as described in [clean.conc.na](#) before working with the BLQ values. The method for handling NA concentrations can affect the output of which points are considered BLQ and which are considered "middle". Values are considered BLQ if they are 0.

`conc.blq` can be set either a scalar indicating what should be done for all BLQ values or a list with elements named "first", "middle", and "last" each set to a scalar.

The meaning of each of the list elements is:

first Values up to the first non-BLQ value. Note that if all values are BLQ, this includes all values.

middle Values that are BLQ between the first and last non-BLQ values.

last Values that are BLQ after the last non-BLQ value

The valid settings for each are:

"drop" Drop the BLQ values

"keep" Keep the BLQ values

a number Set the BLQ values to that number

Value

The concentration and time measurements (data frame) filtered and cleaned as requested relative to BLQ in the middle.

See Also

[clean.conc.na](#)

clean.conc.na	<i>Handle NA values in the concentration measurements as requested by the user.</i>
---------------	-------------------------------------------------------------------------------------

Description

NA concentrations (and their associated times) will be removed then the BLQ values in the middle

Usage

```
clean.conc.na(conc, time, ..., options = list(), conc.na = NULL,  
              check = TRUE)
```

Arguments

conc	Measured concentrations
time	Time of the concentration measurement
...	Additional items to add to the data frame
options	List of changes to the default PKNCA.options for calculations.
conc.na	How to handle NA concentrations? Either 'drop' or a number to impute.
check	Run check.conc.time ?

Value

The concentration and time measurements (data frame) filtered and cleaned as requested relative to NA in the concentration.

exclude	<i>Exclude data points or results from calculations or summarization.</i>
---------	---------------------------------------------------------------------------

Description

Exclude data points or results from calculations or summarization.

Usage

```
exclude(object, reason, mask, FUN)
```

```
## Default S3 method:
```

```
exclude(object, reason, mask, FUN)
```

Arguments

object	The object to exclude data from.
reason	The reason to add as a reason for exclusion.
mask	A logical vector or numeric index of values to exclude (see details).
FUN	A function to operate on the data (one group at a time) to select reasons for exclusions (see details).

Details

Only one of mask or FUN may be given. If FUN is given, it will be called with two arguments: a data.frame (or similar object) that consists of a single group of the data and the full object (e.g. the PKNCAconc object), FUN(current_group, object), and it must return a logical vector equivalent to mask or a character vector with the reason text given when data should be excluded or NA_character_ when the data should be included (for the current exclusion test).

Value

The object with updated information in the exclude column. The exclude column will contain the reason if mask or FUN indicate. If a previous reason for exclusion was given, then subsequent reasons for exclusion will be added to the first with a semicolon space ("; ") separator.

Methods (by class)

- default: The general case for data exclusion

Examples

```
myconc <- PKNCAconc(data.frame(subject=1,
                               time=0:6,
                               conc=c(1, 2, 3, 2, 1, 0.5, 0.25)),
                   conc~time|subject)
exclude(myconc,
```

```
reason="Carryover",
mask=c(TRUE, rep(FALSE, 6)))
```

exclude_nca

Exclude NCA parameters based on examining the parameter set.

Description

Exclude NCA parameters based on examining the parameter set.

Usage

```
exclude_nca_span.ratio(min.span.ratio)

exclude_nca_max.aucinf.pext(max.aucinf.pext)

exclude_nca_min.hl.r.squared(min.hl.r.squared)
```

Arguments

`min.span.ratio` The minimum acceptable span ratio (uses `PKNCA.options("min.span.ratio")` if not provided).

`max.aucinf.pext` The maximum acceptable percent AUC extrapolation (uses `PKNCA.options("max.aucinf.pext")` if not provided).

`min.hl.r.squared` The minimum acceptable r-squared value for half-life (uses `PKNCA.options("min.hl.r.squared")` if not provided).

Functions

- `exclude_nca_span.ratio`: Exclude based on span.ratio
- `exclude_nca_max.aucinf.pext`: Exclude based on AUC percent extrapolated (both observed and predicted)
- `exclude_nca_min.hl.r.squared`: Exclude based on half-life r-squared

Examples

```
my_conc <- PKNCAconc(data.frame(conc=1.1^(3:0),
                                time=0:3,
                                subject=1),
                    conc~time|subject)
my_data <- PKNCAdata(my_conc,
                    intervals=data.frame(start=0, end=Inf,
                                          aucinf.obs=TRUE,
                                          aucpext.obs=TRUE))
my_result <- pk.nca(my_data)
```



```
my_result_excluded <- exclude(my_result,  
                              FUN=exclude_nca_max.aucinf.pext())  
as.data.frame(my_result_excluded)
```

find.tau

Find the repeating interval within a vector of doses

Description

This is intended to find the interval over which x repeats by the rule $\text{unique}(\text{mod}(x, \text{interval}))$ is minimized.

Usage

```
find.tau(x, na.action = na.omit, options = list(), tau.choices = NULL)
```

Arguments

<code>x</code>	the vector to find the interval within
<code>na.action</code>	What to do with NAs in <code>x</code>
<code>options</code>	List of changes to the default PKNCA.options for calculations.
<code>tau.choices</code>	the intervals to look for if the doses are not all equally spaced.

Value

A scalar indicating the repeating interval with the most repetition.

1. If all values are NA then NA is returned.
2. If all values are the same, then 0 is returned.
3. If all values are equally spaced, then that spacing is returned.
4. If one of the choices can minimize the number of unique values, then that is returned.
5. If none of the choices can minimize the number of unique values, then -1 is returned.

findOperator	<i>Find the first occurrence of an operator in a formula and return the left, right, or both sides of the operator.</i>
--------------	-------------------------------------------------------------------------------------------------------------------------

Description

Find the first occurrence of an operator in a formula and return the left, right, or both sides of the operator.

Usage

```
findOperator(x, op, side)
```

Arguments

x	The formula to parse
op	The operator to search for (e.g. +, -, *, /, ...)
side	Which side of the operator would you like to see: 'left', 'right', or 'both'.

Value

The side of the operator requested, NA if requesting the left side of a unary operator, and NULL if the operator is not found.

fit_half_life	<i>Perform the half-life fit given the data. The function simply fits the data without any validation. No selection of points or any other components are done.</i>
---------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------

Description

Perform the half-life fit given the data. The function simply fits the data without any validation. No selection of points or any other components are done.

Usage

```
fit_half_life(data, tlast)
```

Arguments

data	The data to fit. Must have two columns named "log_conc" and "time"
tlast	The time of last observed concentration above the limit of quantification.

Value

A data.frame with one row and columns named "r.squared", "adj.r.squared", "PROB", "lambda.z", "clast.pred", "lambda.z.n.points", "half.life", "span.ratio"

See Also

[pk.calc.half.life](#)

formula.parseFormula *Convert the parsed formula back into the original*

Description

Convert the parsed formula back into the original

Usage

```
## S3 method for class 'parseFormula'
formula(x, drop.groups = FALSE, drop.lhs = FALSE,
  ...)
```

Arguments

x	The parsed formula object to revert to the original
drop.groups	logical. Should the returned formula drop the groups?
drop.lhs	logical. Should the returned formula be one-sided dropping the left hand side?
...	Unused.

Value

A formula (optionally with portions removed)

formula.PKNCAconc *Extract the formula from a PKNCAconc object.*

Description

Extract the formula from a PKNCAconc object.

Usage

```
## S3 method for class 'PKNCAconc'
formula(x, ...)
```

```
## S3 method for class 'PKNCAdose'
formula(x, ...)
```

Arguments

`x` The object to extract the formula from.
`...` Unused

Value

A formula object

<code>formulaOrNames</code>	<i>When parsing the input for PKNCAconc and PKNCAdose allow either formula or character inputs.</i>
-----------------------------	-----------------------------------------------------------------------------------------------------

Description

When parsing the input for PKNCAconc and PKNCAdose allow either formula or character inputs.

Usage

```
formulaOrNames(formula = NA_character_, lhs = NA_character_,
  rhs = NA_character_, groups = c())
```

Arguments

`formula` A formula with an optional left side, right side, and groups after a pipe.
`lhs, rhs, groups` Alternate definition of the formula

Value

A list with named elements for the "lhs", "rhs", and "groups"

<code>geomean</code>	<i>Compute the geometric mean, sd, and CV</i>
----------------------	-----------------------------------------------

Description

Compute the geometric mean, sd, and CV

Usage

```
geomean(x, na.rm = FALSE)
geosd(x, na.rm = FALSE)
geocv(x, na.rm = FALSE)
```

Arguments

x A vector to compute the geometric mean of
na.rm Should missing values be removed?

Value

The scalar value of the geometric mean, geometric standard deviation, or geometric coefficient of variation.

Functions

- `geosd`: Compute the geometric standard deviation.
- `geocv`: Compute the geometric coefficient of variation.

References

Kirkwood T. B.L. Geometric means and measures of dispersion. *Biometrics* 1979; 35: 908-909

`get.best.model` *Extract the best model from a list of models using AIC.list.*

Description

Extract the best model from a list of models using `AIC.list`.

Usage

```
get.best.model(object, ...)
```

Arguments

object the list of models
... Parameters passed to `AIC.list`

Value

The model which is assessed as best. If more than one are equal, the first is chosen.

See Also

[AIC.list](#)

get.first.model *Get the first model from a list of models*

Description

Get the first model from a list of models

Usage

```
get.first.model(object)
```

Arguments

object the list of (lists of, ...) models

Value

The first item in the object that is not a list or NA. If NA is passed in or the list (of lists) is all NA, then NA is returned.

get.interval.cols *Get the columns that can be used in an interval specification*

Description

Get the columns that can be used in an interval specification

Usage

```
get.interval.cols()
```

Value

A list with named elements for each parameter. Each list element contains the parameter definition.

Examples

```
get.interval.cols()
```

get.parameter.deps *Get all columns that depend on a parameter*

Description

Get all columns that depend on a parameter

Usage

```
get.parameter.deps(x)
```

Arguments

x The parameter name (as a character string)

Value

A character vector of parameter names that depend on the parameter x. If none depend on x, then the result will be an empty vector.

getAttributeColumn *Retrieve the value of an attribute column.*

Description

Retrieve the value of an attribute column.

Usage

```
getAttributeColumn(object, attr_name, warn_missing = c("attr", "column"))
```

Arguments

object The object to extract the attribute value from.
attr_name The name of the attribute to extract
warn_missing Give a warning if the "attr"ibute or "column" is missing. Character vector with zero, one, or both of "attr" and "column".

Value

The value of the attribute (or NULL if the attribute is not set or the column does not exist)

getColumnValueOrNot	<i>Get the value from a column in a data frame if the value is a column there, otherwise, the value should be a scalar or the length of the data.</i>
---------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------

Description

Get the value from a column in a data frame if the value is a column there, otherwise, the value should be a scalar or the length of the data.

Usage

```
getColumnValueOrNot(data, value, prefix = "X")
```

Arguments

data	A data.frame or similar object
value	A character string giving the name of a column in the data, a scalar, or a vector the same length as the data
prefix	The prefix to use if a column must be added (it will be used as the full column name if it is not already in the dataset or it will be prepended to the maximum column name if not.)

Value

A list with elements named "data", "name" giving the data with a column named "name" with the value in that column.

getData.PKNCAconc	<i>Extract all the original data from a PKNCAconc or PKNCAdose object</i>
-------------------	---------------------------------------------------------------------------

Description

Extract all the original data from a PKNCAconc or PKNCAdose object

Usage

```
## S3 method for class 'PKNCAconc'
getData(object)
```

```
## S3 method for class 'PKNCAdose'
getData(object)
```

Arguments

object	R object to extract the data from.
--------	------------------------------------

getData.PKNCAdata	<i>Extract all the original data from a PKNCAconc or PKNCAdose object</i>
-------------------	---------------------------------------------------------------------------

Description

Extract all the original data from a PKNCAconc or PKNCAdose object

Usage

```
## S3 method for class 'PKNCAdata'  
getData(object)
```

Arguments

object R object to extract the data from.

getData.PKNCResults	<i>Extract all the original data from a PKNCAconc or PKNCAdose object</i>
---------------------	---------------------------------------------------------------------------

Description

Extract all the original data from a PKNCAconc or PKNCAdose object

Usage

```
## S3 method for class 'PKNCResults'  
getData(object)
```

Arguments

object R object to extract the data from.

getDataName.PKNCAconc *Get the name of the element containing the data for the current object.*

Description

Get the name of the element containing the data for the current object.

Usage

```
## S3 method for class 'PKNCAconc'
getDataName(object)

## S3 method for class 'PKNCAdata'
getDataName(object)

## S3 method for class 'PKNCAdose'
getDataName(object)

## S3 method for class 'PKNCAresults'
getDataName(object)

getDataName(object)

## Default S3 method:
getDataName(object)
```

Arguments

object The object to get the data name from.

Value

A character scalar with the name of the data object (or NULL if the method does not apply).

Methods (by class)

- default: If no data name exists, returns NULL.

getDepVar *Get the dependent variable (left hand side of the formula) from a PKNCA object.*

Description

Get the dependent variable (left hand side of the formula) from a PKNCA object.

Usage

```
getDepVar(x, ...)
```

Arguments

```
x          The object to extract the formula from
...        Unused
```

Value

The vector of the dependent variable from the object.

```
getGroups.PKNCAconc  Get the groups (right hand side after the | from a PKNCA object).
```

Description

Get the groups (right hand side after the | from a PKNCA object).

Usage

```
## S3 method for class 'PKNCAconc'
getGroups(object, form = formula(object), level,
  data = getData(object), sep)

## S3 method for class 'PKNCAdose'
getGroups(...)

## S3 method for class 'PKNCAresults'
getGroups(object, form = formula(object$data$conc),
  level, data = object$result, sep)
```

Arguments

```
object      The object to extract the data from
form        The formula to extract the data from (defaults to the formula from object)
level       optional. If included, this specifies the level(s) of the groups to include. If a
            numeric scalar, include the first level number of groups. If a numeric vector,
            include each of the groups specified by the number. If a character vector, include
            the named group levels.
data        The data to extract the groups from (defaults to the data from object)
sep         Unused (kept for compatibility with the nlme package)
...         Arguments passed to other getGroups functions
```

Value

A data frame with the (selected) group columns.

getIndepVar	<i>Get the independent variable (right hand side of the formula) from a PKNCA object.</i>
-------------	-------------------------------------------------------------------------------------------

Description

Get the independent variable (right hand side of the formula) from a PKNCA object.

Usage

```
getIndepVar(x, ...)
```

Arguments

x	The object to extract the formula from
...	Unused

Value

The vector of the independent variable from the object.

interp.extrap.conc	<i>Interpolate concentrations between measurements or extrapolate concentrations after the last measurement.</i>
--------------------	------------------------------------------------------------------------------------------------------------------

Description

interpolate.conc and extrapolate.conc returns an interpolated (or extrapolated) concentration. interp.extrap.conc will choose whether interpolation or extrapolation is required and will also operate on many concentrations. These will typically be used to estimate the concentration between two measured concentrations or after the last measured concentration. Of note, these functions will not extrapolate prior to the first point.

Usage

```
interp.extrap.conc(conc, time, time.out, lambda.z = NA,
  clast = pk.calc.clast.obs(conc, time), options = list(),
  interp.method = NULL, extrap.method = "AUCinf", ..., conc.blq = NULL,
  conc.na = NULL, check = TRUE)
```

```
interpolate.conc(conc, time, time.out, options = list(),
  interp.method = NULL, conc.blq = NULL, conc.na = NULL,
  conc.origin = 0, ..., check = TRUE)
```

```
extrapolate.conc(conc, time, time.out, lambda.z = NA,
```

```
clast = pk.calc.clast.obs(conc, time), extrap.method = "AUCinf",
options = list(), conc.na = NULL, conc.blq = NULL, ..., check = TRUE)
```

```
interp.extrap.conc.dose(conc, time, time.dose, route.dose = "extravascular",
duration.dose = NA, time.out, out.after = FALSE, options = list(),
conc.blq = NULL, conc.na = NULL, ..., check = TRUE)
```

Arguments

conc	Measured concentrations
time	Time of the concentration measurement
time.out	Time when interpolation is requested (vector for interp.extrap.conc, scalar otherwise)
lambda.z	The elimination rate constant. NA will prevent extrapolation.
clast	The last observed concentration above the limit of quantification. If not given, clast is calculated from pk.calc.clast.obs
options	List of changes to the default PKNCA.options for calculations.
interp.method	The method for interpolation (either 'lin up/log down' or 'linear')
extrap.method	The method for extrapolation: "AUCinf", "AUClast", or "AUCall". See details for usage.
...	Additional arguments passed to interpolate.conc or extrapolate.conc .
conc.blq	How to handle BLQ values. (See clean.conc.blq for usage instructions.)
conc.na	How to handle NA concentrations. (See clean.conc.na)
check	Run check.conc.time , clean.conc.blq , and clean.conc.na ?
conc.origin	The concentration before the first measurement. conc.origin is typically used to set predose values to zero (default), set a predose concentration for endogenous compounds, or set predose concentrations to NA if otherwise unknown.
time.dose	Time of the dose
route.dose	What is the route of administration ("intravascular" or "extravascular"). See the details below for how this parameter is used.
duration.dose	What is the duration of administration? See the details below for how this parameter is used.
out.after	Should interpolation occur from the data before (FALSE) or after (TRUE) the interpolated point? See the details below for how this parameter is used. It only has a meaningful effect at the instant of an IV bolus dose.

Details

extrap.method 'AUCinf' Use lambda.z to extrapolate beyond the last point with the half-life.

'AUCall' If the last point is above the limit of quantification or missing, this is identical to 'AUCinf'. If the last point is below the limit of quantification, then linear interpolation between the Clast and the next BLQ is used for that interval and all additional points are extrapolated as 0.

'AUClast' Extrapolates all points after the last above the limit of quantification as 0.

duration.dose and direction.out are ignored if route.dose == "extravascular". direction.out is ignored if duration.dose > 0.

route.dose and duration.dose affect how interpolation/extrapolation of the concentration occurs at the time of dosing. If route.dose == "intravascular" and duration.dose == 0 then extrapolation occurs for an IV bolus using `pk.calc.c0` with the data after dosing. Otherwise (either route.dose == "extravascular" or duration.dose > 0), extrapolation occurs using the concentrations before dosing and estimating the half-life (or more precisely, estimating λ_z). Finally, direction.out can change the direction of interpolation in cases with route.dose == "intravascular" and duration.dose == 0. When direction.out == "before" interpolation occurs only with data before the dose (as is the case for route.dose == "extravascular"), but if direction.out == "after" interpolation occurs from the data after dosing.

Value

The interpolated or extrapolated concentration value as a scalar float.

Functions

- `interpolate.conc`: Interpolate concentrations through Tlast (inclusive)
- `extrapolate.conc`: Extrapolate concentrations after Tlast
- `interp.extrap.conc.dose`: Interpolate and extrapolate concentrations without interpolating or extrapolating beyond doses.

See Also

[pk.calc.clast.obs](#), [pk.calc.half.life](#), [pk.calc.c0](#)

<code>merge.splitlist</code>	<i>Merge two or more lists with a data.frame 'groupid' attribute defining the matching.</i>
------------------------------	---------------------------------------------------------------------------------------------

Description

Merge two or more lists with a data.frame 'groupid' attribute defining the matching.

Usage

```
## S3 method for class 'splitlist'
merge(...)
```

Arguments

... lists with 'groupid' attributes

Details

The merge is equivalent to a `full_join` where items missing from one or the other item will be missing, but the element(s) will exist.

Value

A list of lists with elements for the matching items between the 'groupid's of each of the input lists. The output list will have a new 'groupid' attribute added with additional columns to indicate the of each input to its output location. If the inputs are named, then each list item will be named the same as the input name.

`model.frame.PKNCAconc` *Extract the columns used in the formula (in order) from a PKNCAconc or PKNCAdose object.*

Description

Extract the columns used in the formula (in order) from a PKNCAconc or PKNCAdose object.

Usage

```
## S3 method for class 'PKNCAconc'  
model.frame(formula, ...)  
  
## S3 method for class 'PKNCAdose'  
model.frame(formula, ...)
```

Arguments

<code>formula</code>	The object to use (parameter name is <code>formula</code> to use the generic function)
<code>...</code>	Unused

Value

A data frame with the columns from the object in formula order.

normalize_exclude	<i>Normalize the exclude column by setting blanks to NA</i>
-------------------	-------------------------------------------------------------

Description

Normalize the exclude column by setting blanks to NA

Usage

```
normalize_exclude(object)
```

Arguments

object The object to extract the exclude column from

Value

The exclude vector where NA indicates not to exclude and anything else indicates to exclude.

parseFormula	<i>Parse a formula into its component parts.</i>
--------------	--------------------------------------------------

Description

This function supports parsing

Usage

```
parseFormula(form, require.groups = FALSE, require.two.sided = FALSE)
```

Arguments

form the formula to extract into its parts
 require.groups is it an error not to have groups?
 require.two.sided is it an error to have a one-sided formula?

Details

This function extracts the left hand side (lhs), right hand side (rhs), groups (groups and as a formula, grpFormula), the environment (env, and the original left/right hand side of the model (model).

This function borrows heavily from the parseGroupFormula function in the doBy package.

Value

A parseFormula class list with elements of

model The left~right side of the model (excluding groups)

lhs The call for the left hand side

rhs The call for the right hand side (excluding groups)

groups The call for the groups

groupFormula A formula form of the groups

env The original formula's environment

Examples

```
parseFormula("a~b", require.groups=FALSE)
## parseFormula("a~b", require.groups=TRUE) # This is an error
parseFormula("a~b|c")
parseFormula("a~b|c")$groups
```

pk.business	<i>Run any function with a maximum missing fraction of X and 0s possibly counting as missing. The maximum fraction missing comes from PKNCA.options("max.missing").</i>
-------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Description

Note that all missing values are removed prior to calling the function. The function is called with the

Usage

```
pk.business(FUN, zero.missing = FALSE, max.missing)
```

Arguments

FUN	function to run. The function is called as FUN(x, ...) with missing values removed.
zero.missing	Are zeros counted as missing? If TRUE then include them in the missing count.
max.missing	The maximum fraction of the data allowed to be missing (a number between 0 and 1, inclusive).

Value

A version of FUN that can be called with parameters that are checked for missingness (and zeros) with missing (and zeros) removed before the call. If max.missing is exceeded, then NA is returned.

pk.calc.ae *Calculate amount excreted (typically in urine or feces)*

Description

Calculate amount excreted (typically in urine or feces)

Usage

```
pk.calc.ae(conc, volume, check = TRUE)
```

Arguments

conc	The concentration in the sample
volume	The volume (or mass) of the sample
check	Should the concentration and volume data be checked?

Details

The units for the concentration and volume should match such that $\text{sum}(\text{conc} \times \text{volume})$ has units of mass or moles.

Value

The amount excreted during the interval

See Also

[pk.calc.clr](#), [pk.calc.fe](#)

pk.calc.aucint *Calculate the AUC over an interval with interpolation and/or extrapolation of concentrations for the beginning and end of the interval.*

Description

Calculate the AUC over an interval with interpolation and/or extrapolation of concentrations for the beginning and end of the interval.

Usage

```
pk.calc.aucint(conc, time, interval = NULL, start = NULL, end = NULL,
  clast = pk.calc.clast.obs(conc, time), lambda.z = NA, time.dose = NULL,
  route = "extravascular", duration.dose = 0, method = NULL,
  auc.type = "AUClast", ..., options = list())
```

```
pk.calc.aucint.last(conc, time, start = NULL, end = NULL, time.dose, ...,
  options = list())
```

```
pk.calc.aucint.all(conc, time, start = NULL, end = NULL, time.dose, ...,
  options = list())
```

```
pk.calc.aucint.inf.obs(conc, time, start = NULL, end = NULL, time.dose,
  lambda.z, clast.obs, ..., options = list())
```

```
pk.calc.aucint.inf.pred(conc, time, start = NULL, end = NULL, time.dose,
  lambda.z, clast.pred, ..., options = list())
```

Arguments

conc	Concentration measured
time	Time of concentration measurement (must be monotonically increasing and the same length as the concentration data)
interval	Numeric vector of two numbers for the start and end time of integration
start, end	The start and end of the interval (cannot be given if interval is given)
clast, clast.obs, clast.pred	The last concentration above the limit of quantification; this is used for AUCinf calculations. If provided as clast.obs (observed clast value, default), AUCinf is AUCinf.obs. If provided as clast.pred, AUCinf is AUCinf.pred.
lambda.z	The elimination rate (in units of inverse time) for extrapolation
time.dose, route, duration.dose	The time of doses, route of administration, and duration of dose used with interpolation and extrapolation of concentration data (see interp.extrap.conc.dose). If NULL, interp.extrap.conc will be used instead (assuming that no doses affecting concentrations are in the interval).
method	The method for integration (either 'lin up/log down' or 'linear')
auc.type	The type of AUC to compute. Choices are 'AUCinf', 'AUClast', and 'AUCall'.
...	Additional arguments passed to <code>pk.calc.auxc</code> and <code>interp.extrap.conc</code>
options	List of changes to the default PKNCA.options for calculations.

Functions

- `pk.calc.aucint.last`: Interpolate or extrapolate concentrations for AUClast
- `pk.calc.aucint.all`: Interpolate or extrapolate concentrations for AUCall
- `pk.calc.aucint.inf.obs`: Interpolate or extrapolate concentrations for AUCinf.obs
- `pk.calc.aucint.inf.pred`: Interpolate or extrapolate concentrations for AUCinf.pred

See Also

[pk.calc.auxc](#), [PKNCA.options](#), [interp.extrap.conc.dose](#)

pk.calc.aucpext *Calculate the AUC percent extrapolated*

Description

Calculate the AUC percent extrapolated

Usage

```
pk.calc.aucpext(auclast, aucinf)
```

Arguments

auclast the area under the curve from time 0 to the last measurement above the limit of quantification

aucinf the area under the curve from time 0 to infinity

Value

The numeric value of the AUC percent extrapolated or NA_real_ if any of the following are true is.na(aucinf), is.na(auclast), aucinf <= 0, or auclast <= 0.

pk.calc.auxc *A compute the Area Under the (Moment) Curve*

Description

Compute the area under the curve (AUC) and the area under the moment curve (AUMC) for pharmacokinetic (PK) data. AUC and AUMC are used for many purposes when analyzing PK in drug development.

Usage

```
pk.calc.auxc(conc, time, interval = c(0, Inf),
  clast = pk.calc.clast.obs(conc, time, check = FALSE), lambda.z = NA,
  auc.type = "AUClast", options = list(), method = NULL,
  conc.blq = NULL, conc.na = NULL, check = TRUE, fun.linear, fun.log,
  fun.inf)
```

```
pk.calc.auc(conc, time, ..., options = list())
```

```
pk.calc.auc.last(conc, time, ..., options = list())
```

```

pk.calc.auc.inf(conc, time, ..., options = list(), lambda.z)
pk.calc.auc.inf.obs(conc, time, clast.obs, ..., options = list(), lambda.z)
pk.calc.auc.inf.pred(conc, time, clast.pred, ..., options = list(), lambda.z)
pk.calc.auc.all(conc, time, ..., options = list())
pk.calc.aumc(conc, time, ..., options = list())
pk.calc.aumc.last(conc, time, ..., options = list())
pk.calc.aumc.inf(conc, time, ..., options = list(), lambda.z)
pk.calc.aumc.inf.obs(conc, time, clast.obs, ..., options = list(), lambda.z)
pk.calc.aumc.inf.pred(conc, time, clast.pred, ..., options = list(), lambda.z)
pk.calc.aumc.all(conc, time, ..., options = list())

```

Arguments

conc	Concentration measured
time	Time of concentration measurement (must be monotonically increasing and the same length as the concentration data)
interval	Numeric vector of two numbers for the start and end time of integration
clast, clast.obs, clast.pred	The last concentration above the limit of quantification; this is used for AUCinf calculations. If provided as clast.obs (observed clast value, default), AUCinf is AUCinf,obs. If provided as clast.pred, AUCinf is AUCinf,pred.
lambda.z	The elimination rate (in units of inverse time) for extrapolation
auc.type	The type of AUC to compute. Choices are 'AUCinf', 'AUClast', and 'AUCall'.
options	List of changes to the default PKNCA.options for calculations.
method	The method for integration (either 'lin up/log down' or 'linear')
conc.blq	How to handle BLQ values in between the first and last above LOQ concentrations. (See clean.conc.blq for usage instructions.)
conc.na	How to handle missing concentration values. (See clean.conc.na for usage instructions.)
check	Run check.conc.time , clean.conc.blq , and clean.conc.na ?
fun.linear	The function to use for integration of the linear part of the curve (not required for AUC or AUMC functions)
fun.log	The function to use for integration of the logarithmic part of the curve (if log integration is used; not required for AUC or AUMC functions)

fun.inf The function to use for extrapolation from the final measurement to infinite time
 (not required for AUC or AUMC functions).
 ... For functions other than pk.calc.auxc, these values are passed to pk.calc.auxc

Details

pk.calc.auc.last is simply a shortcut setting the interval parameter to $c(0, "last")$.
 Extrapolation beyond Clast occurs using the half-life and Clast,obs; Clast,pred is not yet supported.

Value

A numeric value for the AU(M)C

Functions

- pk.calc.auc: Compute the area under the curve
- pk.calc.auc.last: Compute the AUClast.
- pk.calc.auc.inf: Compute the AUCinf
- pk.calc.auc.inf.obs: Compute the AUCinf with the observed Clast.
- pk.calc.auc.inf.pred: Compute the AUCinf with the predicted Clast.
- pk.calc.auc.all: Compute the AUCall.
- pk.calc.aumc: Compute the area under the moment curve
- pk.calc.aumc.last: Compute the AUMClast.
- pk.calc.aumc.inf: Compute the AUMCinf
- pk.calc.aumc.inf.obs: Compute the AUMCinf with the observed Clast.
- pk.calc.aumc.inf.pred: Compute the AUMCinf with the predicted Clast.
- pk.calc.aumc.all: Compute the AUMCall.

References

Gabrielsson J, Weiner D. "Section 2.8.1 Computation methods - Linear trapezoidal rule." Pharmacokinetic & Pharmacodynamic Data Analysis: Concepts and Applications, 4th Edition. Stockholm, Sweden: Swedish Pharmaceutical Press, 2000. 162-4.

Gabrielsson J, Weiner D. "Section 2.8.3 Computation methods - Log-linear trapezoidal rule." Pharmacokinetic & Pharmacodynamic Data Analysis: Concepts and Applications, 4th Edition. Stockholm, Sweden: Swedish Pharmaceutical Press, 2000. 164-7.

See Also

[pk.calc.auc.all](#), [pk.calc.auc.last](#), [clean.conc.blq](#)

Examples

```
myconc <- c(0, 1, 2, 1, 0.5, 0.25, 0)
mytime <- c(0, 1, 2, 3, 4, 5, 6)
pk.calc.auc(myconc, mytime, interval=c(0, 6))
pk.calc.auc(myconc, mytime, interval=c(0, Inf))
```

`pk.calc.c0`*Estimate the concentration at dosing time for an IV bolus dose.*

Description

Estimate the concentration at dosing time for an IV bolus dose.

Usage

```
pk.calc.c0(conc, time, time.dose = 0, method = c("c0", "logslope", "c1"),  
           check = TRUE)
```

```
pk.calc.c0.method.logslope(conc, time, time.dose = 0, check = TRUE)
```

```
pk.calc.c0.method.c0(conc, time, time.dose = 0, check = TRUE)
```

```
pk.calc.c0.method.c1(conc, time, time.dose = 0, check = TRUE)
```

```
pk.calc.c0.method.set0(conc, time, time.dose = 0, check = TRUE)
```

```
pk.calc.c0.method.cmin(conc, time, time.dose = 0, check = TRUE)
```

Arguments

<code>conc</code>	The observed concentrations
<code>time</code>	The observed times
<code>time.dose</code>	The time when dosing occurred
<code>method</code>	The order of methods to test (see details)
<code>check</code>	Check the <code>conc</code> and <code>time</code> inputs

Details

Methods available for interpolation are below, and each has its own specific function.

`c0` If the observed `conc` at `time.dose` is nonzero, return that. This method should always be used first.

`logslope` Compute the semilog line between the first two measured times, and use that line to extrapolate backward to `time.dose`

`c1` Use the first point after `time.dose`

Value

The estimated concentration at time 0.

Functions

- `pk.calc.c0.method.logslope`: Semilog regress the first and second points after `time.dose`. This method will return NA if the second conc after `time.dose` is 0 or greater than the first.
- `pk.calc.c0.method.c0`: Use $C_0 = \text{conc}[\text{time } \text{time.dose}]$ if it is nonzero.
- `pk.calc.c0.method.c1`: Use $C_0 = C_1$.
- `pk.calc.c0.method.set0`: Use $C_0 = 0$ (typically used for single dose oral and IV infusion)
- `pk.calc.c0.method.cmin`: Use $C_0 = C_{\text{min}}$ (typically used for multiple dose oral and IV infusion)

`pk.calc.cav` *Calculate the average concentration during an interval.*

Description

Calculate the average concentration during an interval.

Usage

```
pk.calc.cav(auclast, start, end)
```

Arguments

<code>auclast</code>	The area under the curve during the interval
<code>start</code>	The starting time of the interval
<code>end</code>	The ending time of the interval

Value

The C_{av} (average concentration during the interval)

`pk.calc.ceoi` *Determine the concentration at the end of infusion*

Description

Determine the concentration at the end of infusion

Usage

```
pk.calc.ceoi(conc, time, duration.dose = NA, check = TRUE)
```


Arguments

conc	Concentration measured
time	Time of concentration measurement
duration.dose	The duration for the dosing administration (typically from IV infusion)
check	Run <code>check.conc.time</code> ?

Value

The concentration at the end of the infusion, NA if duration.dose is NA, or NA if all time != duration.dose

pk.calc.cl	<i>Calculate the (observed oral) clearance</i>
------------	------------------------------------------------

Description

Calculate the (observed oral) clearance

Usage

```
pk.calc.cl(dose, auc)
```

Arguments

dose	the dose administered
auc	The area under the concentration-time curve.

Details

If dose is the same length as the other inputs, then the output will be the same length as all of the inputs; the function assumes that you are calculating for multiple intervals simultaneously. If the inputs other than dose are scalars and dose is a vector, then the function assumes multiple doses were given in a single interval, and the sum of the doses will be used for the calculation.

Value

the numeric value of the total (CL) or observed oral clearance (CL/F)

References

Gabrielsson J, Weiner D. "Section 2.5.1 Derivation of clearance." Pharmacokinetic & Pharmacodynamic Data Analysis: Concepts and Applications, 4th Edition. Stockholm, Sweden: Swedish Pharmaceutical Press, 2000. 86-7.

pk.calc.clast.obs	<i>Determine the last observed concentration above the limit of quantification (LOQ).</i>
-------------------	-------------------------------------------------------------------------------------------

Description

If Tlast is NA (due to no non-missing above LOQ measurements), this will return NA.

Usage

```
pk.calc.clast.obs(conc, time, check = TRUE)
```

Arguments

conc	Concentration measured
time	Time of concentration measurement
check	Run check.conc.time?

Value

The last observed concentration above the LOQ

pk.calc.clr	<i>Calculate renal clearance</i>
-------------	----------------------------------

Description

Calculate renal clearance

Usage

```
pk.calc.clr(ae, auc)
```

Arguments

ae	The amount excreted in urine (as a numeric scalar or vector)
auc	The area under the curve (as a numeric scalar or vector)

Details

The units for the ae and auc should match such that ae/auc has units of volume/time.

Value

The renal clearance as a number

See Also

[pk.calc.ae](#), [pk.calc.fe](#)

pk.calc.cmax	<i>Determine maximum observed PK concentration</i>
--------------	----------------------------------------------------

Description

Determine maximum observed PK concentration

Usage

```
pk.calc.cmax(conc, check = TRUE)
```

```
pk.calc.cmin(conc, check = TRUE)
```

Arguments

conc	Concentration measured
check	Run check.conc.time?

Value

a number for the maximum concentration or NA if all concentrations are missing

Functions

- `pk.calc.cmin`: Determine the minimum observed PK concentration

pk.calc.ctrough	<i>Determine the trough (predose) concentration</i>
-----------------	-----------------------------------------------------

Description

Determine the trough (predose) concentration

Usage

```
pk.calc.ctrough(conc, time, start)
```

Arguments

conc	Observed concentrations during the interval
time	Times of conc observations
start	Starting time of the interval

Value

The concentration when time == start. If none match, then NA

pk.calc.deg.fluc *Determine the degree of fluctuation*

Description

Determine the degree of fluctuation

Usage

pk.calc.deg.fluc(cmax, cmin, cav)

Arguments

cmax	The maximum observed concentration
cmin	The minimum observed concentration
cav	The average concentration in the interval

Value

The degree of fluctuation around the average concentration.

pk.calc.dn *Determine dose normalized NCA parameter*

Description

Determine dose normalized NCA parameter

Usage

pk.calc.dn(parameter, dose)

Arguments

parameter	Parameter to dose normalize
dose	Dose in units compatible with the area under the curve

Value

a number for dose normalized AUC

Examples

pk.calc.dn(90, 10)

pk.calc.f *Calculate the absolute (or relative) bioavailability*

Description

Calculate the absolute (or relative) bioavailability

Usage

pk.calc.f(dose1, auc1, dose2, auc2)

Arguments

dose1	The dose administered in route or method 1
auc1	The AUC from 0 to infinity or 0 to tau administered in route or method 1
dose2	The dose administered in route or method 2
auc2	The AUC from 0 to infinity or 0 to tau administered in route or method 2

pk.calc.fe *Calculate fraction excreted (typically in urine or feces)*

Description

Calculate fraction excreted (typically in urine or feces)

Usage

pk.calc.fe(ae, dose)

Arguments

ae	The amount excreted (as a numeric scalar or vector)
dose	The dose (as a numeric scalar or vector)

Details

The units for ae and dose should be the same so that ae/dose is a unitless fraction.

Value

The fraction of dose excreted.

See Also

[pk.calc.ae](#), [pk.calc.clr](#)

pk.calc.half.life *Compute the half-life and associated parameters*

Description

Compute the half-life and associated parameters

Usage

```
pk.calc.half.life(conc, time, tmax, tlast, manually.selected.points = FALSE,
  options = list(), min.hl.points = NULL, adj.r.squared.factor = NULL,
  conc.blq = NULL, conc.na = NULL, first.tmax = NULL,
  allow.tmax.in.half.life = NULL, check = TRUE)
```

Arguments

conc	Concentration measured
time	Time of concentration measurement
tmax	Time of maximum concentration (will be calculated and included in the return data frame if not given)
tlast	Time of last concentration above the limit of quantification (will be calculated and included in the return data frame if not given)
manually.selected.points	Have the points being input been manually selected? The impact of setting this to TRUE is that no selection for the best points will be done. When TRUE, this option causes the options of <code>adj.r.squared.factor</code> , <code>min.hl.points</code> , and <code>allow.tmax.in.half.life</code> to be ignored.
options	List of changes to the default <code>PKNCA.options</code> for calculations.
min.hl.points	The minimum number of points that must be included to calculate the half-life
adj.r.squared.factor	The allowance in adjusted r-squared for adding another point.
conc.blq	See clean.conc.blq
conc.na	See clean.conc.na
first.tmax	See pk.calc.tmax .
allow.tmax.in.half.life	Allow the concentration point for tmax to be included in the half-life slope calculation.
check	Run check.conc.time , clean.conc.blq , and clean.conc.na ?

Details

If `manually.selected.points` is FALSE (default), the half-life is calculated by computing the best fit line for all available sets of points. The best one is chosen by the following rules in order:

- At least `min.hl.points` points included
- A `lambda.z > 0`
- The best adjusted r-squared (within `adj.r.squared.factor`)
- The one with the most points included

If `manually.selected.points` is TRUE, the `conc` and `time` data are used as-is without any form of selection for the best-fit half-life.

Value

A data frame with one row and columns for

tmax Time of maximum observed concentration (only included if not given as an input)

tlast Time of last observed concentration above the LOQ (only included if not given as an input)

r.squared coefficient of determination

adj.r.squared adjusted coefficient of determination

lambda.z elimination rate

lambda.z.time.first first time for half-life calculation

lambda.z.n.points number of points in half-life calculation

clast.pred Concentration at `tlast` as predicted by the half-life line

half.life half-life

span.ratio span ratio [ratio of half-life to time used for half-life calculation]

References

Gabrielsson J, Weiner D. "Section 2.8.4 Strategies for estimation of lambda-z." Pharmacokinetic & Pharmacodynamic Data Analysis: Concepts and Applications, 4th Edition. Stockholm, Sweden: Swedish Pharmaceutical Press, 2000. 167-9.

pk.calc.kel

Calculate the elimination rate (Kel)

Description

Calculate the elimination rate (Kel)

Usage

pk.calc.kel(mrt)

Arguments

mrt the mean residence time

Value

the numeric value of the elimination rate

pk.calc.mrt	<i>Calculate the mean residence time (MRT) for single-dose data or linear multiple-dose data.</i>
-------------	---------------------------------------------------------------------------------------------------

Description

Calculate the mean residence time (MRT) for single-dose data or linear multiple-dose data.

Usage

```
pk.calc.mrt(auc, aumc)
```

```
pk.calc.mrt.iv(auc, aumc, duration.dose)
```

Arguments

auc the AUC from 0 to infinity or 0 to tau

aumc the AUMC from 0 to infinity or 0 to tau

duration.dose The duration of the dose (usually an infusion duration for an IV infusion)

Value

the numeric value of the mean residence time

Functions

- `pk.calc.mrt.iv`: MRT for an IV infusion

See Also

[pk.calc.mrt.md](#)

pk.calc.mrt.md	<i>Calculate the mean residence time (MRT) for multiple-dose data with nonlinear kinetics.</i>
----------------	------------------------------------------------------------------------------------------------

Description

Calculate the mean residence time (MRT) for multiple-dose data with nonlinear kinetics.

Usage

```
pk.calc.mrt.md(auctau, aumctau, aucinf, tau)
```

Arguments

auctau	the AUC from time 0 to the end of the dosing interval (tau).
aumctau	the AUMC from time 0 to the end of the dosing interval (tau).
aucinf	the AUC from time 0 to infinity (typically using single-dose data)
tau	the dosing interval

Details

Note that if `aucinf == auctau` (as would be the assumption with linear kinetics), the equation becomes the same as the single-dose MRT.

See Also

[pk.calc.mrt](#)

pk.calc.ptr	<i>Determine the peak-to-trough ratio</i>
-------------	-------------------------------------------

Description

Determine the peak-to-trough ratio

Usage

```
pk.calc.ptr(cmax, cmin)
```

Arguments

cmax	The maximum observed concentration
cmin	The minimum observed concentration

Value

The ratio of `cmax` to `cmin` (if `cmin == 0`, NA)

`pk.calc.swing`*Determine the PK swing*

Description

Determine the PK swing

Usage

```
pk.calc.swing(cmax, cmin)
```

Arguments

<code>cmax</code>	The maximum observed concentration
<code>cmin</code>	The minimum observed concentration

Value

The swing above the minimum concentration. If `cmin` is zero, then the result is infinity.

`pk.calc.thalf.eff`*Calculate the effective half-life*

Description

Calculate the effective half-life

Usage

```
pk.calc.thalf.eff(mrt)
```

Arguments

<code>mrt</code>	the mean residence time to infinity
------------------	-------------------------------------

Value

the numeric value of the effective half-life

pk.calc.tlag	<i>Determine the observed lag time (time before the first concentration above the limit of quantification or above the first concentration in the interval)</i>
--------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------

Description

Determine the observed lag time (time before the first concentration above the limit of quantification or above the first concentration in the interval)

Usage

```
pk.calc.tlag(conc, time)
```

Arguments

conc	The observed concentrations
time	The observed times

Value

The time associated with the first increasing concentration

pk.calc.tlast	<i>Determine time of last observed concentration above the limit of quantification.</i>
---------------	-----------------------------------------------------------------------------------------

Description

NA will be returned if all conc are NA or 0.

Usage

```
pk.calc.tlast(conc, time, check = TRUE)
```

```
pk.calc.tfirst(conc, time, check = TRUE)
```

Arguments

conc	Concentration measured
time	Time of concentration measurement
check	Run <code>check.conc.time?</code>

Value

The time of the last observed concentration measurement

Functions

- `pk.calc.tfirst`: Determine the first concentration above the limit of quantification.

`pk.calc.tmax`*Determine time of maximum observed PK concentration*

Description

Input restrictions are:

1. the `conc` and `time` must be the same length,
2. the `time` may have no NAs,

NA will be returned if:

1. the length of `conc` and `time` is 0
2. all `conc` is 0 or NA

Usage

```
pk.calc.tmax(conc, time, options = list(), first.tmax = NULL,  
             check = TRUE)
```

Arguments

<code>conc</code>	Concentration measured
<code>time</code>	Time of concentration measurement
<code>options</code>	List of changes to the default PKNCA.options for calculations.
<code>first.tmax</code>	If there is more than time that matches the maximum concentration, should the first be considered as Tmax? If not, then the last is considered Tmax.
<code>check</code>	Run check.conc.time ?

Value

the time of the maximum concentration

pk.calc.vd	<i>Calculate the volume of distribution (Vd) or observed volume of distribution (Vd/F)</i>
------------	--------------------------------------------------------------------------------------------

Description

Calculate the volume of distribution (Vd) or observed volume of distribution (Vd/F)

Usage

```
pk.calc.vd(dose, aucinf, lambda.z)
```

Arguments

dose	One or more doses given during an interval
aucinf	Area under the curve to infinity (either predicted or observed).
lambda.z	Elimination rate constant

Details

If dose is the same length as the other inputs, then the output will be the same length as all of the inputs; the function assumes that you are calculating for multiple intervals simultaneously. If the inputs other than dose are scalars and dose is a vector, then the function assumes multiple doses were given in a single interval, and the sum of the doses will be used for the calculation.

Value

The observed volume of distribution

pk.calc.vss	<i>Calculate the steady-state volume of distribution (Vss)</i>
-------------	----------------------------------------------------------------

Description

Calculate the steady-state volume of distribution (Vss)

Usage

```
pk.calc.vss(cl, mrt)
```

Arguments

cl	the clearance
mrt	the mean residence time

Value

the volume of distribution at steady-state

pk.calc.vz	<i>Calculate the terminal volume of distribution (Vz)</i>
------------	-----------------------------------------------------------

Description

Calculate the terminal volume of distribution (Vz)

Usage

```
pk.calc.vz(cl, lambda.z)
```

Arguments

cl	the clearance (or apparent observed clearance)
lambda.z	the elimination rate

pk.nca	<i>Compute NCA parameters for each interval for each subject.</i>
--------	-------------------------------------------------------------------

Description

The computations assume that all calculation options are set in [PKNCA.options](#).

Usage

```
pk.nca(data)
```

Arguments

data	A PKNCAdata object
------	--------------------

Value

A data frame with a row for each interval for each grouping in the concentration data.

See Also

[PKNCAdata](#), [PKNCA.options](#)

pk.nca.interval *Compute all PK parameters for a single concentration-time data set*

Description

For one subject/time range, compute all available PK parameters. All the internal options should be set by `PKNCA.options` prior to running. The only part that changes with a call to this function is the concentration and time.

Usage

```
pk.nca.interval(conc, time, volume, duration.conc, dose, time.dose,
  duration.dose, route, conc.group = NULL, time.group = NULL,
  volume.group = NULL, duration.conc.group = NULL, dose.group = NULL,
  time.dose.group = NULL, duration.dose.group = NULL, route.group = NULL,
  include_half.life = NULL, exclude_half.life = NULL, interval,
  options = list())
```

Arguments

`conc`, `conc.group`
Concentration measured for the current interval or all data for the group

`time`, `time.group`
Time of concentration measurement for the current interval or all data for the group

`volume`, `volume.group`
The volume (or mass) of the concentration measurement for the current interval or all data for the group (typically for urine and fecal measurements)

`duration.conc`, `duration.conc.group`
The duration of the concentration measurement for the current interval or all data for the group (typically for urine and fecal measurements)

`dose`, `dose.group`
Dose amount (may be a scalar or vector) for the current interval or all data for the group

`time.dose`, `time.dose.group`
Time of the dose for the current interval or all data for the group (must be the same length as dose or dose.group)

`duration.dose`, `duration.dose.group`
The duration of the dose administration for the current interval or all data for the group (typically zero for extravascular and intravascular bolus and nonzero for intravascular infusion)

`route`, `route.group`
The route of dosing for the current interval or all data for the group

`include_half.life`
An optional boolean vector of the concentration measurements to include in the half-life calculation. If given, no half-life point selection will occur.

exclude_half.life	An optional boolean vector of the concentration measurements to exclude from the half-life calculation.
interval	One row of an interval definition (see check.interval.specification for how to define the interval).
options	List of changes to the default PKNCA.options for calculations.

Value

A data frame with the start and end time along with all PK parameters for the interval

See Also

[check.interval.specification](#)

pk.tss	<i>Compute the time to steady-state (tss)</i>
--------	-----------------------------------------------

Description

Compute the time to steady-state (tss)

Usage

```
pk.tss(..., type = c("monoexponential", "stepwise.linear"), check = TRUE)
```

Arguments

...	Passed to pk.tss.monoexponential or pk.tss.stepwise.linear .
type	The type of Tss to calculate, either stepwise.linear or monoexponential
check	See pk.tss.data.prep

Value

A data frame with columns as defined from [pk.tss.monoexponential](#) and/or [pk.tss.stepwise.linear](#).

See Also

[pk.tss.monoexponential](#), [pk.tss.stepwise.linear](#)

pk.tss.data.prep *Clean up the time to steady-state parameters and return a data frame for use by the tss calculators.*

Description

Clean up the time to steady-state parameters and return a data frame for use by the tss calculators.

Usage

```
pk.tss.data.prep(conc, time, subject, treatment, subject.dosing, time.dosing,
  options = list(), conc.blq = NULL, conc.na = NULL, check = TRUE, ...)
```

Arguments

conc	Concentration measured
time	Time of concentration measurement
subject	Subject identifiers (used as a random effect in the model)
treatment	Treatment description (if missing, all subjects are assumed to be on the same treatment)
subject.dosing	Subject number for dosing
time.dosing	Time of dosing
options	List of changes to the default <code>PKNCA.options</code> for calculations.
conc.blq	See clean.conc.blq
conc.na	See clean.conc.na
check	Run check.conc.time?
...	Discarded inputs to allow generic calls between tss methods.

Value

a data frame with columns for concentration, time, subject, and treatment.

pk.tss.monoexponential *Compute the time to steady state using nonlinear, mixed-effects modeling of trough concentrations.*

Description

Trough concentrations are selected as concentrations at the time of dosing. An exponential curve is then fit through the data with a different magnitude by treatment (as a factor) and a random steady-state concentration and time to steady-state by subject (see `random.effects` argument).

Usage

```
pk.tss.monoexponential(..., tss.fraction = 0.9, output = c("population",
  "popind", "individual", "single"), check = TRUE, verbose = FALSE)
```

Arguments

...	See pk.tss.data.prep
tss.fraction	The fraction of steady-state required for calling steady-state
output	Which types of outputs should be produced? <code>population</code> is the population estimate for time to steady-state (from an nlme model), <code>popind</code> is the individual estimate (from an nlme model), <code>individual</code> fits each individual separately with a gnls model (requires more than one individual; use <code>single</code> for one individual), and <code>single</code> fits all the data to a single gnls model.
check	See pk.tss.data.prep .
verbose	Describe models as they are run, show convergence of the model (passed to the nlme function), and additional details while running.

Value

A scalar float for the first time when steady-state is achieved or NA if it is not observed.

References

Maganti L, Panebianco DL, Maes AL. Evaluation of Methods for Estimating Time to Steady State with Examples from Phase 1 Studies. *AAPS Journal* 10(1):141-7. doi:10.1208/s12248-008-9014-y

See Also

[pk.tss](#), [pk.tss.stepwise.linear](#)

pk.tss.monoexponential.individual

A helper function to estimate individual and single outputs for mono-exponential time to steady-state.

Description

This function is not intended to be called directly. Please use `pk.tss.monoexponential`.

Usage

```
pk.tss.monoexponential.individual(data, output = c("individual", "single"),
  verbose = FALSE)
```

Arguments

data	a data frame as prepared by pk.tss.data.prep . It must contain at least columns for subject, time, conc, and tss.constant.
output	a character vector requesting the output types.
verbose	Show verbose output.

Details

If no model converges, then the `tss.monoexponential.single` and/or `tss.monoexponential.individual` column will be set to NA.

Value

A data frame with either one row (if population output is provided) or one row per subject (if popind is provided). The columns will be named `tss.monoexponential.population` and/or `tss.monoexponential.popind`.

`pk.tss.monoexponential.population`

A helper function to estimate population and popind outputs for monoexponential time to steady-state.

Description

This function is not intended to be called directly. Please use `pk.tss.monoexponential`.

Usage

```
pk.tss.monoexponential.population(data, output = c("population", "popind"),
  verbose = FALSE)
```

Arguments

data	a data frame as prepared by pk.tss.data.prep . It must contain at least columns for subject, time, conc, and tss.constant.
output	a character vector requesting the output types.
verbose	Show verbose output.

Details

If no model converges, then the `tss.monoexponential.population` column will be set to NA. If the best model does not include a random effect for subject on Tss then the `tss.monoexponential.popind` column of the output will be set to NA.

Value

A data frame with either one row (if population output is provided) or one row per subject (if popind is provided). The columns will be named `tss.monoexponential.population` and/or `tss.monoexponential.popind`.

`pk.tss.stepwise.linear`

Compute the time to steady state using stepwise test of linear trend

Description

A linear slope is fit through the data to find when it becomes non-significant. Note that this is less preferred than the `pk.tss.monoexponential` due to the fact that with more time or more subjects the performance of the test changes (see reference).

Usage

```
pk.tss.stepwise.linear(..., min.points = 3, level = 0.95, verbose = FALSE,  
  check = TRUE)
```

Arguments

<code>...</code>	See pk.tss.data.prep
<code>min.points</code>	The minimum number of points required for the fit
<code>level</code>	The confidence level required for assessment of steady-state
<code>verbose</code>	Describe models as they are run, show convergence of the model (passed to the <code>nlme</code> function), and additional details while running.
<code>check</code>	See pk.tss.data.prep

Details

The model is fit with a different magnitude by treatment (as a factor, if given) and a random slope by subject (if given). A minimum of `min.points` is required to fit the model.

Value

A scalar float for the first time when steady-state is achieved or NA if it is not observed.

References

Maganti L, Panebianco DL, Maes AL. Evaluation of Methods for Estimating Time to Steady State with Examples from Phase 1 Studies. *AAPS Journal* 10(1):141-7. doi:10.1208/s12248-008-9014-y

See Also

[pk.tss.monoexponential](#)

Description

Compute pharmacokinetic (PK) noncompartmental analysis (NCA) parameters.

Details

A common workflow would load data from a file or database into a data.frame then run the following

Examples

```
## Not run:
# Load concentration-time data into a data.frame called d.conc
# with columns named "conc", "time", and "subject".
my.conc <- PKNCAconc(d.conc, conc~time|subject)
# Load dose-time data into a data.frame called d.dose
# with columns named "dose", "time", and "subject".
my.dose <- PKNCAdose(d.dose, dose~time|subject)
# Combine the concentration-time and dose-time data into an object
# ready for calculations.
my.data <- PKNCAdata(my.conc, my.dose)
# Perform the calculations
my.results <- pk.nca(my.data)
# Look at summary results
summary(my.results)
# Look at a listing of results
as.data.frame(my.results)

## End(Not run)
```

PKNCA.choose.option *Choose either the value from an option list or the current set value for an option.*

Description

Choose either the value from an option list or the current set value for an option.

Usage

```
PKNCA.choose.option(name, value = NULL, options = list())
```

Arguments

name	The option name requested.
value	A value to check for the option (NULL to choose not to check the value).
options	The non-default options to choose from.

Value

The value of the option first from the `options` list and if it is not there then from the current settings.

See Also

[PKNCA.options](#)

PKNCA.options	<i>Set default options for PKNCA functions</i>
---------------	------------------------------------------------

Description

This function will set the default PKNCA options. If given no inputs, it will provide the current option set. If given name/value pairs, it will set the option (as in the [options](#) function). If given a name, it will return the value for the parameter. If given the default option as true, it will provide the default options.

Usage

```
PKNCA.options(..., default = FALSE, check = FALSE, name, value)
```

Arguments

...	options to set or get the value for
default	(re)sets all default options
check	check a single option given, but do not set it (for validation of the values when used in another function)
name	An option name to use with the value.
value	An option value (paired with the name) to set or check (if NULL,).

Details

Options are either for calculation or summary functions. Calculation options are required for a calculation function to report a result (otherwise the reported value will be NA). Summary options are used during summarization and are used for assessing what values are included in the summary.

See the vignette 'Options for Controlling PKNCA' for a current list of options.

Value

If...

no arguments are given returns the current options.

a value is set (including the defaults) returns NULL

a single value is requested the current value of that option is returned as a scalar

multiple values are requested the current values of those options are returned as a list

See Also

[PKNCA.choose.option](#)

Examples

```
PKNCA.options()  
PKNCA.options(default=TRUE)  
PKNCA.options("auc.method")  
PKNCA.options(name="auc.method")  
PKNCA.options(auc.method="lin up/log down", min.hl.points=3)
```

PKNCA.options.describe

Describe a PKNCA.options option by name.

Description

Describe a PKNCA.options option by name.

Usage

```
PKNCA.options.describe(name)
```

Arguments

name The option name requested.

Value

A character string of the description.

See Also

[PKNCA.options](#)

PKNCA.set.summary *Define how NCA parameters are summarized.*

Description

Define how NCA parameters are summarized.

Usage

```
PKNCA.set.summary(name, point, spread, rounding = list(signif = 3),  
  reset = FALSE)
```

Arguments

name	The parameter name or a vector of parameter names. It must have already been defined (see add.interval.col).
point	The function to calculate the point estimate for the summary. The function will be called as <code>point(x)</code> and must return a scalar value (typically a number, NA, or a string).
spread	Optional. The function to calculate the spread (or variability). The function will be called as <code>spread(x)</code> and must return a scalar or two-long vector (typically a number, NA, or a string).
rounding	Instructions for how to round the value of point and spread. It may either be a list or a function. If it is a list, then it must have a single entry with a name of either "signif" or "round" and a value of the digits to round. If a function, it is expected to return a scalar number or character string with the correct results for an input of either a scalar or a two-long vector.
reset	Reset all the summary instructions

Value

All current summary settings (invisibly)

See Also

[summary.PKNCAresults](#)

 PKNCAconc

 Create a PKNCAconc object

Description

Create a PKNCAconc object

Usage

```
PKNCAconc(data, ...)

## Default S3 method:
PKNCAconc(data, ...)

## S3 method for class 'tbl_df'
PKNCAconc(data, ...)

## S3 method for class 'data.frame'
PKNCAconc(data, formula, subject, time.nominal, exclude,
  duration, volume, exclude_half.life, include_half.life, ...)
```

Arguments

data	A data frame with concentration (or amount for urine/feces), time, and the groups defined in formula.
...	Ignored.
formula	The formula defining the concentration~time groups or amount~time groups for urine/feces (In the remainder of the documentation, "concentration" will be used to describe concentration or amount.)
subject	The column indicating the subject number (used for plotting). If not provided, this defaults to the beginning of the inner groups: For example with concentration~time Study+Subject the inner groups start with the first grouping variable before a /, Subject. If there is only one grouping variable, it is assumed to be the subject (e.g. concentration~time Subject), and if there are multiple grouping variables without a /, subject is assumed to be the last one. For single-subject data, it is assigned as NULL.
time.nominal	(optional) The name of the nominal time column (if the main time variable is actual time. The time.nominal is not used during calculations; it is available to assist with data summary and checking.
exclude	(optional) The name of a column with concentrations to exclude from calculations and summarization. If given, the column should have values of NA or "" for concentrations to include and non-empty text for concentrations to exclude.
duration	(optional) The duration of collection as is typically used for concentration measurements in urine or feces.

volume (optional) The volume (or mass) of collection as is typically used for urine or feces measurements.

exclude_half.life, include_half.life Points to exclude from the half-life calculation (still using normal selection rules for the other points) or to include for the half-life (using specifically those points and bypassing automatic point selection).

Value

A PKNCAconc object that can be used for automated NCA.

See Also

[PKNCAdata](#), [PKNCAdose](#)

PKNCAdata

Create a PKNCAdata object.

Description

PKNCAdata combines PKNCAconc and PKNCAdose and adds in the intervals for PK calculations.

Usage

```
PKNCAdata(data.conc, data.dose, ...)

## S3 method for class 'PKNCAconc'
PKNCAdata(data.conc, data.dose, ...)

## S3 method for class 'PKNCAdose'
PKNCAdata(data.conc, data.dose, ...)

## Default S3 method:
PKNCAdata(data.conc, data.dose, ..., formula.conc,
           formula.dose, intervals, options = list())
```

Arguments

data.conc Concentration data as a PKNCAconc object or a data frame

data.dose Dosing data as a PKNCAdose object (see details)

... arguments passed to PKNCAdata.default

formula.conc Formula for making a PKNCAconc object with data.conc. This must be given if data.conc is a data.frame, and it must not be given if data.conc is a PKNCAconc object.

formula.dose Formula for making a PKNCAdose object with data.dose. This must be given if data.dose is a data.frame, and it must not be given if data.dose is a PKNCAdose object.

intervals	A data frame with the AUC interval specifications as defined in check.interval.specification . If missing, this will be automatically chosen by choose.auc.intervals . (see details)
options	List of changes to the default PKNCA.options for calculations.

Details

If `data.dose` is not given or is NA, then the `intervals` must be given. At least one of `data.dose` and `intervals` must be given.

Value

A PKNCAdata object with concentration, dose, interval, and calculation options stored (note that PKNCAdata objects can also have results after a NCA calculations are done to the data).

See Also

[PKNCAconc](#), [PKNCAdose](#), [choose.auc.intervals](#), [pk.nca](#)

 PKNCAdose

Create a PKNCAdose object

Description

Create a PKNCAdose object

Usage

```
PKNCAdose(data, ...)

## Default S3 method:
PKNCAdose(data, ...)

## S3 method for class 'tbl_df'
PKNCAdose(data, ...)

## S3 method for class 'data.frame'
PKNCAdose(data, formula, route, rate, duration,
           time.nominal, exclude, ...)
```

Arguments

data	A data frame with time and the groups defined in formula.
...	Ignored.
formula	The formula defining the <code>dose.amount~time groups</code> where <code>time</code> is the time of the dosing and <code>dose.amount</code> is the amount administered at that time (see Details).

route	Define the route of administration. The value may be either a column name from the data (checked first) or a character string of either "extravascular" or "intravascular" (checked second). If given as a column name, then every value of the column must be either "extravascular" or "intravascular".
rate, duration	(optional) for "intravascular" dosing, the rate or duration of dosing. If given as a character string, it is the name of a column from the data, and if given as a number, it is the value for all doses. Only one may be given, and if neither is given, then the dose is assumed to be a bolus (duration=0). If rate is given, then the dose amount must be given (the left hand side of the formula).
time.nominal	(optional) The name of the nominal time column (if the main time variable is actual time. The time.nominal is not used during calculations; it is available to assist with data summary and checking.
exclude	(optional) The name of a column with concentrations to exclude from calculations and summarization. If given, the column should have values of NA or "" for concentrations to include and non-empty text for concentrations to exclude.

Details

The formula for a PKNCAdose object can be given three ways: one-sided (missing left side), one-sided (missing right side), or two-sided. Each of the three ways can be given with or without groups. When given one-sided missing the left side, the left side can either be omitted or can be given as a period (.): `~time|treatment+subject` and `.~time|treatment+subject` are identical, and dose-related NCA parameters will all be reported as not calculable (for example, clearance). When given one-sided missing the right side, the right side must be specified as a period (.): `dose~.|treatment+subject`, and only a single row may be given per group. When the right side is missing, PKNCA assumes that the same dose is given in every interval. When given as a two-sided formula

Value

A PKNCAconc object that can be used for automated NCA.

See Also

[PKNCAconc](#), [PKNCAdata](#)

PKNCAresults

Generate a PKNCAresults object

Description

This function should not be run directly. The object is created for summarization and plotting.

Usage

```
PKNCAresults(result, data, exclude)
```

Arguments

result	a data frame with NCA calculation results and groups. Each row is one interval and each column is a group name or the name of an NCA parameter.
data	The PKNCAdata used to generate the result
exclude	(optional) The name of a column with concentrations to exclude from calculations and summarization. If given, the column should have values of NA or "" for concentrations to include and non-empty text for concentrations to exclude.

Value

A PKNCAresults object with each of the above within.

plot.PKNCAconc	<i>Plot a PKNCAconc object</i>
----------------	--------------------------------

Description

Plot a PKNCAconc object

Usage

```
## S3 method for class 'PKNCAconc'
plot(x, ..., groups = x$subject,
      panel.formula = parseFormula(x)$groupFormula, panel.formula.update)

## S3 method for class 'PKNCAdata'
plot(x, ...)
```

Arguments

x	The object to plot
...	Additional arguments passed to xyplot
groups	The grouping variable for the plot (typically the subject column)
panel.formula	The formula used for the call to xyplot (defaults to the group formula of x)
panel.formula.update	Updates to the panel.formula to simplify modifications without having to fully specify the formula.

Value

A trellis object of the plot(s)

```
print.PKNCAconc      Print and/or summarize a PKNCAconc or PKNCAdose object.
```

Description

Print and/or summarize a PKNCAconc or PKNCAdose object.

Usage

```
## S3 method for class 'PKNCAconc'
print(x, n = 6, summarize = FALSE, ...)

## S3 method for class 'PKNCAconc'
summary(object, n = 0, summarize = TRUE, ...)

## S3 method for class 'PKNCAdose'
print(x, n = 6, summarize = FALSE, ...)

## S3 method for class 'PKNCAdose'
summary(object, n = 0, summarize = TRUE, ...)
```

Arguments

x	The object to print
n	The number of rows of data to show (see head)
summarize	Summarize the nested number of groups
...	Arguments passed to <code>print.formula</code> and <code>print.data.frame</code>
object	The object to summarize

```
print.PKNCAdata      Print a PKNCAdata object
```

Description

Print a PKNCAdata object

Usage

```
## S3 method for class 'PKNCAdata'
print(x, ...)
```

Arguments

x	The object to print
...	Arguments passed on to print.PKNCAconc and print.PKNCAdose

print.provenance	<i>Print the summary of a provenance object</i>
------------------	-------------------------------------------------

Description

Print the summary of a provenance object

Usage

```
## S3 method for class 'provenance'
print(x, ...)
```

Arguments

x	The object to be printed
...	Ignored

Value

invisible text of the printed information

roundingSummarize	<i>During the summarization of PKNCAresults, do the rounding of values based on the instructions given.</i>
-------------------	-------------------------------------------------------------------------------------------------------------

Description

During the summarization of PKNCAresults, do the rounding of values based on the instructions given.

Usage

```
roundingSummarize(x, name)
```

Arguments

x	The values to summarize
name	The NCA parameter name (matching a parameter name in PKNCA.set.summary)

Value

A string of the rounded value

roundString	<i>Round a value to a defined number of digits printing out trailing zeros, if applicable.</i>
-------------	------------------------------------------------------------------------------------------------

Description

Round a value to a defined number of digits printing out trailing zeros, if applicable.

Usage

```
roundString(x, digits = 0, si_range = Inf)
```

Arguments

x	The number to round
digits	integer indicating the number of decimal places
si_range	See help for signifString (and you likely want to round with signifString if you want to use this argument)

Details

Values that are not standard numbers like Inf, NA, and NaN are returned as "Inf", "NA", and NaN.

Value

A string with the value

See Also

[round](#), [signifString](#)

setAttributeColumn	<i>Add an attribute to an object where the attribute is added as a name to the names of the object.</i>
--------------------	---------------------------------------------------------------------------------------------------------

Description

Add an attribute to an object where the attribute is added as a name to the names of the object.

Usage

```
setAttributeColumn(object, attr_name, col_or_value, col_name, default_value,
  stop_if_default, warn_if_default, message_if_default)
```


Arguments

object	The object to set the attribute column on.
attr_name	The attribute name to set
col_or_value	If this exists as a column in the data, it is used as the col_name. If not, this becomes the default_value.
col_name	The name of the column within the dataset to use (if missing, uses attr_name)
default_value	The value to fill in the column if the column does not exist (the column is filled with NA if it does not exist and no value is provided).
stop_if_default, warn_if_default, message_if_default	A character string to provide as an error, a warning, or a message to the user if the default_value is used. They are tested in order (if stop, the code stops; if warning, the message is ignored; and message last).

Value

The object with the attribute column added to the data.

See Also

[getAttributeColumn](#)

setDuration	<i>Set the duration of dosing or measurement</i>
-------------	--------------------------------------------------

Description

Set the duration of dosing or measurement

Usage

```
setDuration(object, ...)

## S3 method for class 'PKNCAdose'
setDuration(object, duration, rate, dose, ...)
```

Arguments

object	An object to set a duration on
...	Arguments passed to another setDuration function
duration	The value to set for the duration or the name of the column in the data to use for the duration.
rate	(for PKNCAdose objects only) The rate of infusion
dose	(for PKNCAdose objects only) The dose amount

Value

The object with duration set

setExcludeColumn	<i>Set the exclude parameter on an object</i>
------------------	-----------------------------------------------

Description

This function adds the exclude column to an object. To change the exclude value, use the [exclude](#) function.

Usage

```
setExcludeColumn(object, exclude, dataname = "data")
```

Arguments

object	The object to set the exclude column on.
exclude	The column name to set as the exclude value.
dataname	The name of the data.frame within the object to add the exclude column to.

Value

The object with an exclude column and attribute

setRoute	<i>Set the dosing route</i>
----------	-----------------------------

Description

Set the dosing route

Usage

```
setRoute(object, ...)

## S3 method for class 'PKNCAdose'
setRoute(object, route, ...)
```

Arguments

object	A PKNCAdose object
...	Arguments passed to another setRoute function
route	A character string indicating one of the following: the column from the data which indicates the route of administration, a scalar indicating the route of administration for all subjects, or a vector indicating the route of administration for each dose in the dataset.

Value

The object with an updated route

signifString	<i>Round a value to a defined number of significant digits printing out trailing zeros, if applicable.</i>
--------------	------------------------------------------------------------------------------------------------------------

Description

Round a value to a defined number of significant digits printing out trailing zeros, if applicable.

Usage

```
signifString(x, digits = 6, si_range = 6)
```

```
## S3 method for class 'data.frame'  
signifString(x, digits = 6, si_range = 6)
```

```
## Default S3 method:  
signifString(x, digits = 6, si_range = 6)
```

Arguments

x	The number to round
digits	integer indicating the number of significant digits
si_range	integer (or Inf) indicating when to switch to scientific notation instead of floating point. Zero indicates always use scientific; Inf indicates to never use scientific notation; otherwise, scientific notation is used when $\text{abs}(\log_{10}(x)) > \text{si_range}$.

Details

Values that are not standard numbers like Inf, NA, and NaN are returned as "Inf", "NA", and NaN.

Value

A string with the value

See Also

[signif](#), [roundString](#)

```
sort.interval.cols      Sort the interval columns by dependencies.
```

Description

Columns are always to the right of columns that they depend on.

Usage

```
## S3 method for class 'interval.cols'
sort()
```

```
split.PKNCAconc      Divide into groups
```

Description

split.PKNCAconc divides data into individual groups defined by [getGroups.PKNCAconc](#).

Usage

```
## S3 method for class 'PKNCAconc'
split(x, f = getGroups(x), drop = TRUE, ...)

## S3 method for class 'PKNCAdata'
split(x, ...)

## S3 method for class 'PKNCAdose'
split(x, f = getGroups(x), drop = TRUE, ...)
```

Arguments

x	the object to split
f	the groups to use for splitting the object
drop	logical indicating if levels that do not occur should be dropped.
...	Ignored.

Details

If x is NA then a list with NA as the only element and a "groupid" attribute of an empty data.frame is returned.

Value

A list of objects with an attribute of groupid consisting of a data.frame with columns for each group.

summary.PKNCAdata	<i>Summarize a PKNCAdata object showing important details about the concentration, dosing, and interval information.</i>
-------------------	--------------------------------------------------------------------------------------------------------------------------

Description

Summarize a PKNCAdata object showing important details about the concentration, dosing, and interval information.

Usage

```
## S3 method for class 'PKNCAdata'
summary(object, ...)
```

Arguments

object	The PKNCAdata object to summarize.
...	arguments passed on to print.PKNCAdata

summary.PKNCAresults	<i>Summarize PKNCA results</i>
----------------------	--------------------------------

Description

Summarize PKNCA results

Usage

```
## S3 method for class 'PKNCAresults'
summary(object, ...,
        drop.group = object$data$conc$subject, summarize.n.per.group = TRUE,
        not.requested.string = ".", not.calculated.string = "NC")
```

Arguments

object	The results to summarize
...	Ignored.
drop.group	Which group(s) should be dropped from the formula?
summarize.n.per.group	Should a column for N be added (TRUE or FALSE)? Note that N is maximum number of parameter results for any parameter; if no parameters are requested for a group, then N will be NA.
not.requested.string	A character string to use when a parameter summary was not requested for a parameter within an interval.

not.calculated.string

A character string to use when a parameter summary was requested, but the point estimate AND spread calculations (if applicable) returned NA.

Details

Excluded results will not be included in the summary.

Value

A data frame of NCA parameter results summarized according to the summarization settings.

See Also

[PKNCA.set.summary](#)

superposition	<i>Compute noncompartmental superposition for repeated dosing</i>
---------------	-------------------------------------------------------------------

Description

Compute noncompartmental superposition for repeated dosing

Usage

```
superposition(conc, ...)

## S3 method for class 'PKNCAconc'
superposition(conc, ...)

## S3 method for class 'numeric'
superposition(conc, time, dose.input, tau, dose.times = 0,
  dose.amount, n.tau = Inf, options = list(), lambda.z,
  clast.pred = FALSE, tlast, additional.times = c(), check.blq = TRUE,
  interp.method = NULL, extrap.method = "AUCinf",
  steady.state.tol = 0.001, ...)
```

Arguments

conc	Concentration measured
...	Additional arguments passed to the half.life function if required to compute lambda.z.
time	Time of concentration measurement
dose.input	The dose given to generate the conc and time inputs. If missing, output doses will be assumed to be equal to the input dose.
tau	The dosing interval

dose.times	The time of dosing within the dosing interval. The $\min(\text{dose.times})$ must be ≥ 0 , and the $\max(\text{dose.times})$ must be $< \tau$. There may be more than one dose times given as a vector.
dose.amount	The doses given for the output. Linear proportionality will be used from the input to output if they are not equal. The length of dose.amount must be either 1 or matching the length of dose.times.
n.tau	The number of tau dosing intervals to simulate or Inf for steady-state.
options	The PKNCA.options to use for the calculation (passed on to subsequent functions like pk.calc.half.life).
lambda.z	The elimination rate (from the half life calculation, used to extrapolate beyond the maximum time observed).
clast.pred	To use predicted as opposed to observed Clast, either give the value for clast.pred here or set it to true (for automatic calculation from the half-life).
tlast	The time of last observed concentration above the limit of quantification. This is calculated if not provided.
additional.times	Times to include in the final outputs in addition to the standard times (see details). All $\min(\text{additional.times})$ must be ≥ 0 , and the $\max(\text{additional.times})$ must be $\leq \tau$.
check.blq	Must the first concentration measurement be below the limit of quantification?
interp.method	See interp.extrap.conc
extrap.method	See interp.extrap.conc
steady.state.tol	The tolerance for assessing if steady-state has been achieved (between 0 and 1, exclusive).

Details

The returned superposition times will include all of the following times: 0 (zero), dose.times, time modulo tau (shifting time for each dose time as well), additional.times, and tau.

Value

A data frame with columns named "conc" and "time".

See Also

[interp.extrap.conc](#)

```
tss.monoexponential.generate.formula
```

A helper function to generate the formula and starting values for the parameters in monoexponential models.

Description

A helper function to generate the formula and starting values for the parameters in monoexponential models.

Usage

```
tss.monoexponential.generate.formula(data)
```

Arguments

data The data used for the model

Value

a list with elements for each of the variables

Index

add.interval.col, 4, 64
addProvenance, 6, 12
adj.r.squared, 6
AIC.list, 7, 21
as.data.frame, 7
as.data.frame.PKNCAresults, 7

business.cv (business.mean), 8
business.geocv (business.mean), 8
business.geomean (business.mean), 8
business.max (business.mean), 8
business.mean, 8
business.median (business.mean), 8
business.min (business.mean), 8
business.range (business.mean), 8
business.sd (business.mean), 8

check.conc.time, 9, 13, 14, 29, 37, 41–43,
46, 51, 52, 57
check.conversion, 10
check.interval.deps, 10, 11
check.interval.specification, 10, 11, 12,
56, 67
checkProvenance, 6, 11
choose.auc.intervals, 12, 67
clean.conc.blq, 13, 29, 37, 38, 46, 57
clean.conc.na, 13, 14, 14, 29, 37, 46, 57

exclude, 15, 74
exclude_nca, 16
exclude_nca_max.aucinf.pext
(exclude_nca), 16
exclude_nca_min.hl.r.squared
(exclude_nca), 16
exclude_nca_span.ratio (exclude_nca), 16
extrapolate.conc (interp.extrap.conc),
28

find.tau, 12, 17
findOperator, 18

fit_half_life, 18
formula.parseFormula, 19
formula.PKNCAconc, 19
formula.PKNCAdose (formula.PKNCAconc),
19
formulaOrNames, 20

geocv (geomean), 20
geomean, 20
geosd (geomean), 20
geosd, (geomean), 20
get.best.model, 7, 21
get.first.model, 22
get.interval.cols, 11, 22
get.parameter.deps, 11, 23
getAttributeColumn, 23, 73
getColumnValueOrNot, 24
getData.PKNCAconc, 24
getData.PKNCAdata, 25
getData.PKNCAdose (getData.PKNCAconc),
24
getData.PKNCAresults, 25
getDataName (getDataName.PKNCAconc), 26
getDataName.PKNCAconc, 26
getDepVar, 26
getGroups.PKNCAconc, 27, 76
getGroups.PKNCAdose
(getGroups.PKNCAconc), 27
getGroups.PKNCAresults
(getGroups.PKNCAconc), 27
getIndepVar, 28

head, 70

interp.extrap.conc, 28, 35, 79
interp.extrap.conc.dose, 35, 36
interpolate.conc (interp.extrap.conc),
28

merge.splitlist, 30

- model.frame.PKNCAconc, 31
- model.frame.PKNCAdose
 - (model.frame.PKNCAconc), 31
- normalize_exclude, 32
- options, 62
- parseFormula, 32
- pk.business, 33
- pk.calc.ae, 34, 43, 45
- pk.calc.auc, 12
- pk.calc.auc(pk.calc.auxc), 36
- pk.calc.auc.all, 38
- pk.calc.auc.last, 38
- pk.calc.aucint, 34
- pk.calc.aucpext, 36
- pk.calc.aumc, 12
- pk.calc.aumc(pk.calc.auxc), 36
- pk.calc.auxc, 36, 36
- pk.calc.c0, 30, 39
- pk.calc.cav, 40
- pk.calc.ceoi, 40
- pk.calc.cl, 41
- pk.calc.clast.obs, 29, 30, 42
- pk.calc.clr, 34, 42, 45
- pk.calc.cmax, 43
- pk.calc.cmin(pk.calc.cmax), 43
- pk.calc.ctrough, 43
- pk.calc.deg.fluc, 44
- pk.calc.dn, 44
- pk.calc.f, 45
- pk.calc.fe, 34, 43, 45
- pk.calc.half.life, 12, 19, 30, 46
- pk.calc.kel, 47
- pk.calc.mrt, 48, 49
- pk.calc.mrt.md, 48, 49
- pk.calc.ptr, 49
- pk.calc.swing, 50
- pk.calc.tfirfirst(pk.calc.tlast), 51
- pk.calc.thalf.eff, 50
- pk.calc.tlag, 51
- pk.calc.tlast, 51
- pk.calc.tmax, 46, 52
- pk.calc.vd, 53
- pk.calc.vss, 53
- pk.calc.vz, 54
- pk.nca, 54, 67
- pk.nca.interval, 55
- pk.tss, 56, 58
- pk.tss.data.prep, 56, 57, 58–60
- pk.tss.monoexponential, 56, 57, 60
- pk.tss.monoexponential.individual, 58
- pk.tss.monoexponential.population, 59
- pk.tss.stepwise.linear, 56, 58, 60
- PKNCA, 61
- PKNCA-package (PKNCA), 61
- PKNCA.choose.option, 61, 63
- PKNCA.options, 12–14, 17, 29, 35–37, 46, 52, 54–57, 62, 62, 63, 67
- PKNCA.options.describe, 63
- PKNCA.set.summary, 64, 71, 78
- PKNCAconc, 65, 67, 68
- PKNCAdata, 54, 66, 66, 68
- PKNCAdose, 66, 67, 67
- PKNCAresults, 68
- plot.PKNCAconc, 69
- plot.PKNCAdata(plot.PKNCAconc), 69
- print.PKNCAconc, 70, 70
- print.PKNCAdata, 70, 77
- print.PKNCAdose, 70
- print.PKNCAdose(print.PKNCAconc), 70
- print.provenance, 71
- round, 72
- roundingSummarize, 71
- roundString, 72, 75
- setAttributeColumn, 72
- setDuration, 73
- setExcludeColumn, 74
- setRoute, 74
- signif, 75
- signifString, 72, 75
- sort.interval.cols, 76
- split.PKNCAconc, 76
- split.PKNCAdata(split.PKNCAconc), 76
- split.PKNCAdose(split.PKNCAconc), 76
- summary.PKNCAconc(print.PKNCAconc), 70
- summary.PKNCAdata, 77
- summary.PKNCAdose(print.PKNCAconc), 70
- summary.PKNCAresults, 64, 77
- superposition, 78
- tss.monoexponential.generate.formula, 80